



アルペン カーバイドバー

Round

Model No.	1/4	1/2	1	2	3	4	5	6	7	8
DIA. (mm)	0.50	0.70	0.80	1.00	1.20	1.40	1.60	1.80	2.10	2.30

Inverted Cone

Model No.	33 1/2	34	35	36	37
DIA. (mm)	0.70	0.80	1.00	1.20	1.40

Fissure

Model No.	56	57	58
DIA. (mm)	0.90	1.00	1.20

Tapered Fissure

Model No.	169	170	171
DIA. (mm)	0.90	1.00	1.20

Amalgam Prep

Model No.	256
DIA. (mm)	0.80

Pear

Model No.	330	331	332
DIA. (mm)	0.80	1.00	1.20

Cross Cut Fissure

Model No.	556	557	558	559	560
DIA. (mm)	0.90	1.00	1.20	1.40	1.60

Round End Cross Cut Fissure

Model No.	1556	1557	1558
DIA. (mm)	0.90	1.00	1.20

Cross Cut Tapered Fissure

Model No.	699	700	701	702	703
DIA. (mm)	0.90	1.00	1.20	1.60	2.10

Round End Fissure

Model No.	1156	1157	1158
DIA. (mm)	0.90	1.00	1.20

Round End Tapered Fissure

Model No.	1170	1171
DIA. (mm)	1.00	1.20

Long Tapered Fissure

Model No.	169L	170L	171L
DIA. (mm)	0.90	1.00	1.20

Long Cross Cut Tapered Fissure

Model No.	701L
DIA. (mm)	1.20